



# BioBiteCorrector® Sleep Appliance



## The digital Sleep Appliance

- 3D-Design
- high-quality milled splints
- thermoactive flexibility due to Clearsplint material
- MMA-free (0% residual monomer content)
- Hinge made entirely of titanium
- biocompatible and hypoallergenic
- Infinitely variable protrusion or with spacer rings
- very high wearing and speech comfort
- via our BBC laboratory as BEL contractual dental health insurance benefit possible



## The digital laboratory

## orthodontics and sleep medicine

BBC Orthotec GmbH  
Innsbrucker Str. 2  
83435 Bad Reichenhall  
GERMANY

Phone: +49 (0)8651 / 9650099  
Fax: +49 (0)8651 / 9650098  
e-Mail: [info@bbc-orthotec.de](mailto:info@bbc-orthotec.de)  
Url: [bbc-orthotec.de](http://bbc-orthotec.de)



## The digital laboratory



[bbc-orthotec.de](http://bbc-orthotec.de)

## digital Herbst

- 3D-Design
- without welding seams
- with BBC hinge is no construction bite required
- Available from BBC Laboratories as BEL possible as a statutory health insurance benefit



### Choice of 3 hinges



BioBiteCorrector®  
Hinge



Dentaurem  
Hinge



Scheu  
Hinge


## BioBiteCorrector® TWIN

- no construction bite required
- without occlusal lateral bite blocks
- best wearing comfort and speech comfort
- no constriction of the tongue space
- BioBiteCorrector® hinges reusable on multibracket appliance
- Available from BBC Laboratories as BEL possible as a statutory health insurance benefit



## digital Hyrax-RME

- 3D-Design with anatomical fit
- Flat design




## digital Tongue grid

- 3D-Design
- Anatomical design
- individual design of the spurs (round/pointed) possible



## digital Drilling template




## Pasin-Pin® System

**Distalisation and mesialisation appliance** with 1,1 mm round wire



### Hybrid-RME



**RME**  
with 4 Pins



### BioDisplacementSpring

for the adjustment of misaligned teeth

