

Laboratory Order Form


Digitale BBC-Sleep Appliance



Orthodontic & Dental Sleep Medicine Laboratory

Dental practice:

Company name: _____
 Contact person: _____
 Street: _____
 Postal code/City: _____
 E-Mail: _____
 Phone: _____
 Customer No.: _____
 Reference: _____



Client:

Name: _____

Insurance: legally privately

BEMA Nr.: _____

Your documents:

- Plaster model STL-File Impression

A construction bite is required to manufacture an anti-snoring splint!
 If you transmit the data to us digitally, the buccal scan must be carried out in the construction bite position!

Titratibility:

- Default setting:**
 infinitely variable from -1.5 mm retrusion
 to +2.5 mm protrusion
 Additional protrusion with spacers is possible.
- individual stepless adjustment (4 mm range):
 from mm retrusion
 to mm protrusion

Splint thickness:

- Default splint thickness:**
 1,5 mm
- Individual splint thickness:
 UJ mm (from 1,2 mm to 2,0 mm possible)
 LJ mm (from 1,2 mm to 2,0 mm possible)

Incisal skeletonization in the upper jaw:

- Default without Skeletonization** Individually with incisal skeletonization UJ

Optionen:

- Titanium spacer 2,0 mm (10 pieces) + 18,90 €
 BBC Open End Wrench (2 pieces) + 34,00 €
 Intraoral non-latex elastics Ø 5/16" (8 mm), 2,7 oz. (100 pieces) + 2,00 €
 Intraoral non-latex elastics Ø 1/4" (5 mm), 2,7 oz. (100 pieces) + 2,00 €

Production time:

- Standard 10-12 working days
 Express 5 working days + 80,00 € incl. Express shipping costs
 The production period starts from the day of order until 12:00 p.m.
 excl. shipping time (depending on the delivery conditions of the delivery country)

Notice:

The laboratory order is carried out under the general Terms and Conditions of business, delivery and payment of the dental technician trade.
 All prices exclusive of VAT.